


# Step by step Installation and use of NanoFinder

Author: Kamyar Merabai: [kamyarm@ethz.ch](mailto:kamyarm@ethz.ch)

## Installation:

1. Please go to the link below and download the **Setup.exe** and **NanoFinder examples for nanoparticle analysis.zip**

Link: <https://github.com/ethz-tema/TEMAsingleparticle>



TEMAsingleparticle/NanoFinder/ · ethz-tema/TEMAsingleparticle

Code Issues Pull requests Actions Projects Security Insights

master / TEMAsingleparticle / NanoFinder / Source code /

kammk NanoFinder 5.5.3 Setup and sample analysis ... 2 hours ago History

..

Setup NanoFinder 5.5.3 Setup and sample analysis 2 hours ago

NanoFinder.m Update NanoFinder.m 4 days ago

Here you can also **NanoFinder User guide.pdf** for setting up your analysis and insight on processed data.

TEMAsingleparticle/NanoFinder/ · ethz-tema/TEMAsingleparticle

Code Issues Pull requests Actions Projects Security Insights

master / TEMAsingleparticle / NanoFinder / Source code / Setup /


kammk NanoFinder 5.5.3 Setup and sample analysis ... 2 hours ago History

..

NanoFinder 5.5 User Guide.pdf User Guide for NanoFinder 5.5.3 updated 12 hours ago



NanoFinder examples for nanoparti... Test sets 2 hours ago



Setup.exe NanoFinder 5.5.3 15 hours ago





## 2. Run the setup.exe










3. Search for NanoFinder in your programs directory and copy a shortcut of your program in Desktop or any other location.




## Processing the test sample:


Example test data could be obtained from:

<https://polybox.ethz.ch/index.php/s/jjsf2ozlUjEiC6t>


1. Open the NanoFinder program
2. A terminal window will show up shortly after




3. Wait until the program graphical user interface (GUI) open. It might take a few seconds to minutes depending on the processing power of your computer.





4. Fill in the GUI as shown in the NanoFinder test folder PowerPoint file. Here we try to process the **Test 1 HDF5 file and read setting from GUI** as well, so the GUI will look like the following after we fill it correctly




5. After you fill the GUI, press Run. Then a file browser will popup. direct the file browser to the location of your raw data directory, where you also copied the MDEF.xlsx.





6. Here we want to process the **Test1 raw data as HDF5 file**. In that folder go to the **Before** file and open **MDEF.xlsx**





## 7. Program start processing as following screenshots is showing.


The terminal window, indicator lights in GUI, and popup chart will inform you about the latest analysis that is being done.


\*While the program is processing, it is highly recommended to prevent using Excel, since the program is using it and error may be raised.


In case any error happens, see the terminal error and read about the error more carefully. In case the error is not solved, an Email may be sent with **Nanofinder run status** file attached to it. The Nanofinder run status is generated as soon as you open the Nanofinder program and it is saved at the location of your program. It contains the same data as shown in the terminal window.



















8. After all processing is done, you could access the processed data in the same directory as your raw data input.



9. For more information and a guide on how to set up your analysis please go to the GitHub folder and find the Nanofinder user guide in the same directory as Setup.exe.

References:

1. K. Mehrabi, R. Kaegi, D. Gunther and A. Gundlach-Graham, Emerging investigator series: Automated Single-Nanoparticle Quantification and Classification: A Holistic Study of Particles into and out of Wastewater Treatment Plants in Switzerland, *Environ. Sci.: Nano*, 2021, DOI: 10.1039/D0EN01066A.
2. A. Gundlach-Graham, L. Hendriks, K. Mehrabi and D. Gunther, Monte Carlo Simulation of Low-Count Signals in Time-of-Flight Mass Spectrometry and Its Application to Single-Particle Detection, *Anal. Chem.*, 2018, **90**, 11847-11855.
3. A. Gundlach-Graham and K. Mehrabi, Monodisperse microdroplets: a tool that advances single-particle ICP-MS measurements, *J. Anal. At. Spectrom.*, 2020, **35**, 1727-1739.
4. K. Mehrabi, D. Gunther and A. Gundlach-Graham, Single-particle ICP-TOFMS with online microdroplet calibration for the simultaneous quantification of diverse nanoparticles in complex matrices, *Environmental Science-Nano*, 2019, **6**, 3349-3358.